


Symbol	Symbolname	Symbol description
	Manufacturer	Indicates the manufacturer of the medical device.
	Date of Manufacture	Indicates the date the medical device was manufactured.
	Importer	Indicates the company that imports the product.
	Distributor	Indicates the company that distributes the product.
	Authorised Representative in the EU	Indicates the authorised representative in the European Union.
	Medical Device	Indicates that this product is a medical device.
	Sterilised with Ethylene Oxide	Indicates the device was sterilised using ethylene oxide (EO gas).
	Single Sterile Barrier System	Indicates a single sterile barrier system (solid line).

List of used Symbols

Symbol	Symbolname	Symbol description
	Use-by Date	Indicates the date after which the product must not be used.
	Catalogue Number	Indicates the product's catalogue or reference number.
	CE Marking (with Notified Body Number)	Indicates compliance with Directive 93/42/EEC & Regulation 2017/745; 'XXXX' is the notified body number.
	Unique Device Identification (UDI)	Indicates barcodes and numbers for traceability and product identification. (01) EAN (GS1/HIBCC) (10) LOT (11) date of manufacture (17) Use-by Date
	Caution	Indicates caution is required when handling the product.
	Consult (Electronic) Instructions for Use	Indicates to consult the (electronic) instructions for use.

Symbol	Symbolname	Symbol description
	Batch Code (LOT Number)	Indicates the manufacturer's batch code (LOT) for product identification.
	Do Not Reuse	Indicates the product is for single use only and must not be reused.
	Do Not Resterilise	Indicates the product must not be resterilised.
	Do Not Use if Package is Damaged	Indicates not to use the product if the sterile packaging is damaged.
	Protect from Sunlight	Indicates the product must be protected from direct sunlight.
	Keep Dry	Indicates the product must be kept dry and protected from moisture.