

Article number	Description	Cutting / preparation	Grasping / holding	Clamping / hemostasis / compression	Spreading / retraction	Dilatation / speculum	Suturing / needle holding	Curettage / scraping	Probing / irrigation / aspiration	Removal of staples / sutures	Foreign body manipulation (e.g., Magill)	Intended purpose
15M-580H-145P	Single-use perineal scissors nach Braun-Stadler, 14,5 cm	X										Single-use instrument, sterile. For cutting and preparing tissue.
17M-114H-145P	Single-use artery forceps nach Kocher, straight, 1x2 Zähne, 14,5 cm			X								Single-use instrument, sterile. For clamping, compression, and hemostasis of vessels/tissue.
27M-291H-130P	Single-use nasal speculum nach Kilian, 30 mm shaft					X						Single-use instrument, sterile. For dilatation and exposure of anatomical structures.
1417M-003H-120P	Single-use button-ended cannula, straight, Ø1.2x80 mm									X		Single-use instrument, sterile. For probing, irrigating, or applying/aspirating fluids.
1417M-005H-120P	Single-use wire scissors, serrated/Minge, 12 cm	X										Single-use instrument, sterile. For cutting and preparing tissue.
1417M-004H-150P	Single-use button-ended cannula, curved, Ø1.5x80 mm									X		Single-use instrument, sterile. For probing, irrigating, or applying/aspirating fluids.
1417M-004H-120P	Single-use button-ended cannula, curved, Ø1.2x80 mm									X		Single-use instrument, sterile. For probing, irrigating, or applying/aspirating fluids.
1417M-003H-150P	Single-use button-ended cannula, straight, Ø1.5x80 mm									X		Single-use instrument, sterile. For probing, irrigating, or applying/aspirating fluids.
15M-027H-145P	Single-use scissors, surgical, st/st, curved, 14,5 cm	X										Single-use instrument, sterile. For cutting and preparing tissue.
15M-139H-145P	Single-use scissors nach Metzenbaum, fine, st/st, straight, 14,5 cm	X										Single-use instrument, sterile. For cutting and preparing tissue.
18M-001H-180P	Single-use needle holder, Mayo-Hegar, straight, 18 cm						X					Single-use instrument, sterile. For grasping and guiding suturing needles during wound closure.
26M-352H-210P	Single-use sharp spoon curette, combination, nach Volkmann, rundoval, large, 21 cm							X				Single-use instrument, sterile. For curettage/scraping of tissue.
16M-056H-145P	Single-use forceps, surgical, 1x2 Zähne, fine, 14,5 cm		X									Single-use instrument, sterile. For grasping and holding tissue, materials, or instruments.
15M-133H-145P	Single-use scissors nach Metzenbaum, st/st, curved, 14,5 cm	X										Single-use instrument, sterile. For cutting and preparing tissue.
15M-133H-115P	Single-use scissors nach Metzenbaum, fine, st/st, curved, 11,5 cm	X										Single-use instrument, sterile. For cutting and preparing tissue.
18M-001H-160P	Single-use needle holder, Mayo-Hegar, straight, 16 cm						X					Single-use instrument, sterile. For grasping and guiding suturing needles during wound closure.
19M-174H-165P	Single-use wound hook set nach Senn-Miller, 16,5 cm				X							Single-use instrument, sterile. Zum Spreizen und Retraktieren von Gewebe.
19M-731H-130P	Single-use wound retractor (Weitläner) scharf, with ratchet, double-ended serrated, 3x4 Zähne, 13 cm				X							Single-use instrument, sterile. Zum Spreizen und Retraktieren von Gewebe.
11M-025H-160P	Single-use Magill forceps, small, mit Greifbacken, 17 cm		X									XSingle-use instrument, sterile. For grasping and holding tissue, materials, or instruments. For manipulation/extraction of foreign bodies, e.g., in the airway.
11M-025H-180P	Single-use Magill forceps, large, mit Greifbacken, 18 cm		X									XSingle-use instrument, sterile. For grasping and holding tissue, materials, or instruments. For manipulation/extraction of foreign bodies, e.g., in the airway.
17M-006H-125P	Single-use Micro-mosquito clamp, anatomical, straight, 12,5 cm			X								Single-use instrument, sterile. For clamping, compression, and hemostasis of vessels/tissue.
17M-011H-130P	Single-use mosquito clamp, nach Halsted, straight, 13 cm			X								Single-use instrument, sterile. For clamping, compression, and hemostasis of vessels/tissue.
17M-007H-125P	Single-use mosquito clamp, curved, fine, 12,5 cm			X								Single-use instrument, sterile. For clamping, compression, and hemostasis of vessels/tissue.
17M-317H-165P	Single-use tubing clamp, smooth jaws, 16,5 cm (Priesdylerian)			X								Single-use instrument, sterile. For clamping, compression, and hemostasis of vessels/tissue.
17M-075H-160P	Single-use artery forceps, Rochester Pean, anatomical, curved, 16 cm			X								Single-use instrument, sterile. For clamping, compression, and hemostasis of vessels/tissue.
17M-073H-180P	Single-use artery forceps Rochester Pean, anatomical straight, 18 cm			X								Single-use instrument, sterile. For clamping, compression, and hemostasis of vessels/tissue.
16M-409H-090P	Single-use Splitterpinzette (Falkenfeld), 9 cm		X									Single-use instrument, sterile. For grasping and holding tissue, materials, or instruments.
16M-234H-120P	Single-use Micro-forceps, anatomical, straight, Adson, 12 cm		X									Single-use instrument, sterile. For grasping and holding tissue, materials, or instruments.
16M-055H-145P	Single-use forceps, anatomical, straight, fine, 14,5 cm		X									Single-use instrument, sterile. For grasping and holding tissue, materials, or instruments.
15M-589H-160P	Single-use dressing scissors nach Lister, 16 cm	X										Single-use instrument, sterile. For cutting and preparing tissue.
15M-238H-115CP	Single-use iris scissors, sp/sp, curved, 11,5 cm	X										Single-use instrument, sterile. For cutting and preparing tissue.
16M-043H-145P	Single-use scissors, surgical, sp/st, curved, 14,5 cm	X										Single-use instrument, sterile. For cutting and preparing tissue.
16M-018H-300P	Single-use forceps, anatomical, straight, 30 cm		X									Single-use instrument, sterile. For grasping and holding tissue, materials, or instruments.
16M-004H-120P	forceps surgical straight 1x2 Zähne, 12 cm		X									Single-use instrument, sterile. For grasping and holding tissue, materials, or instruments.
16M-055H-115P	Single-use forceps, anatomical, straight, fine 11,5cm		X									Single-use instrument, sterile. For grasping and holding tissue, materials, or instruments.
17M-335H-200P	Single-use sponge forceps (Pean Murphy-Zange), with ratchet, 20 cm		X	X								Single-use instrument, sterile. For grasping and holding tissue, materials, or instruments. For clamping, compression, and hemostasis of vessels/tissue.
15M-037H-160P	Single-use umbilical cord scissors, rosa Markierung, 16 cm	X										Single-use instrument, sterile. For cutting and preparing tissue.
16M-400H-115P	Single-use forceps, anatomical, straight, fine, 11,5 cm		X									Single-use instrument, sterile. For grasping and holding tissue, materials, or instruments.
15M-003H-160P	Single-use forceps nach Jansen bapnett, curved, 16 cm		X									Single-use instrument, sterile. For grasping and holding tissue, materials, or instruments.
16M-403H-145P	Single-use forceps anatomical straight, fine		X									Single-use instrument, sterile. For grasping and holding tissue, materials, or instruments.
16M-005H-130P	Single-use forceps nach Semken, surgical, 1x2 Zähne, 13 cm		X									Single-use instrument, sterile. For grasping and holding tissue, materials, or instruments.
15M-589H-140P	Single-use dressing scissors nach Lister, 14 cm	X										Single-use instrument, sterile. For cutting and preparing tissue.
15M-589H-180P	Single-use dressing scissors nach Lister, 18 cm	X										Single-use instrument, sterile. For cutting and preparing tissue.
18M-001H-125P	Single-use needle holder, Mayo-Hegar, straight, 12,5 cm						X					Single-use instrument, sterile. For grasping and guiding suturing needles during wound closure.
17M-335H-160P	Single-use sponge forceps (Pean Murphy-Zange), with ratchet, 16 cm		X	X								Single-use instrument, sterile. For grasping and holding tissue, materials, or instruments. For clamping, compression, and hemostasis of vessels/tissue.
14M-17H-003P	Single-use staple remover, steril, 10 cm										X	Single-use instrument, sterile. For removing skin staples and/or suture material.
15M-003H-140P	Single-use scissors, surgical, sp/sp, straight, 14 cm	X										Single-use instrument, sterile. For cutting and preparing tissue.
18M-001H-150P	Single-use needle holder, Mayo-Hegar, straight, 15 cm						X					Single-use instrument, sterile. For grasping and guiding suturing needles during wound closure.
16M-050H-120P	Single-use forceps, anatomical, straight, Adson, 12 cm		X									Single-use instrument, sterile. For grasping and holding tissue, materials, or instruments.
M1417H-001P	Single-use button-ended cannula, 15cm, steril									X		Single-use instrument, sterile. For probing, irrigating, or applying/aspirating fluids.
26M-352H-170P	Single-use sharp spoon curette, combination, nach Volkmann, rundoval, large, 17 cm							X				Single-use instrument, sterile. For curettage/scraping of tissue.
26M-351H-160P	Single-use sharp spoon curette, combination, nach Volkmann, rundoval, medium, 16 cm							X				Single-use instrument, sterile. For curettage/scraping of tissue.
26M-350H-160P	Single-use sharp spoon curette, combination nach Volkmann, rundoval, small, 16 cm							X				Single-use instrument, sterile. For curettage/scraping of tissue.
15M-589H-145P	Single-use dressing scissors Lister, 14,5cm	X										Single-use instrument, sterile. For cutting and preparing tissue.
15M-019H-145P	Single-use scissors surgical, straight, st/st, 14,5cm	X										Single-use instrument, sterile. For cutting and preparing tissue.
14M-17H-002P	Single-use forceps, DeBakey, anatomical, straight, 15,5 cm		X									Single-use instrument, sterile. For grasping and holding tissue, materials, or instruments.
16M-236H-120P	Single-use Micro-forceps, surgical, 1x2 Zähne, Adson, 12 cm		X									Single-use instrument, sterile. For grasping and holding tissue, materials, or instruments.
18M-001H-140P	Single-use needle holder, Mayo-Hegar, straight, 14 cm						X					Single-use instrument, sterile. For grasping and guiding suturing needles during wound closure.
16M-040H-140P	Single-use forceps anatomical straight, 14cm		X									Single-use instrument, sterile. For grasping and holding tissue, materials, or instruments.
15M-035H-145P	Single-use scissors, surgical, sp/st, straight, 14,5 cm	X										Single-use instrument, sterile. For cutting and preparing tissue.
16M-040H-145P	Single-use forceps, anatomical, straight, 14,5 cm		X									Single-use instrument, sterile. For grasping and holding tissue, materials, or instruments.
15M-238H-115P	Single-use iris scissors, sp/sp, straight, 11,5 cm	X										Single-use instrument, sterile. For cutting and preparing tissue.
16M-409H-115P	Single-use Splitterpinzette (Falkenfeld), 11,5 cm		X									Single-use instrument, sterile. For grasping and holding tissue, materials, or instruments.
17M-073H-140P	Single-use artery forceps, Rochester Pean, anatomical, straight, 14 cm			X								Single-use instrument, sterile. For clamping, compression, and hemostasis of vessels/tissue.
16M-004H-145P	Single-use forceps, surgical, 1x2 Zähne, 14,5 cm		X									Single-use instrument, sterile. For grasping and holding tissue, materials, or instruments.